




The Perfect PerCryo Probe For Any Procedure

PERC-24


Diameter 2.4mm
Shaft Length 15cm



	DIAMETER (MM)	LENGTH (MM)
ICEBALL	37	56
-20 C ISOTHERM	24	44
-40 C ISOTHERM	16	36

PERC-24L


Diameter 2.4mm
Shaft Length 23cm



	DIAMETER (MM)	LENGTH (MM)
ICEBALL	37	56
-20 C ISOTHERM	24	44
-40 C ISOTHERM	16	36

PERC-17


Diameter 1.7mm
Shaft Length 15cm




	DIAMETER (MM)	LENGTH (MM)
ICEBALL	33	54
-20 C ISOTHERM	21	42
-40 C ISOTHERM	14	35

PERC-15

Diameter 1.7mm
Shaft Length 15cm



	DIAMETER (MM)	LENGTH (MM)
ICEBALL	32	34
-20 C ISOTHERM	18	20
-40 C ISOTHERM	13	15



*Data collected using gelatin formula which approximates performance (±5mm) in soft tissue at 10 minutes.



HealthTronics
Total Urology. Total Care.

(888) 252-6575 healthtronics.com



Endocare is a subsidiary of HealthTronics®